

MSE672–Lab 1: Basic Imaging

John Dunlap and Gerd Duscher

February 10, 2021

1 OBJECTIVES

- Knowledge of system components
- Basic alignment of TEM
- Acquisition of images

2 SYSTEM

- Operational status and vacuum system of TEM
- Identification of electro-optical parts of TEM
- Basic functionality of computer control

3 HANDS ON EXERCISE

1. Focus of sample – gold nanoparticles on carbon
2. Changing the illumination system
3. Aligning *Illumination Tilt*
4. Changing the projector system
5. Acquire images with different magnifications (mark clearly the indicated magnification)
6. Measure (ruler) the distance between middle of objective lens and objective aperture

4 HOMEWORK I

Calibration of magnification and rotation of a set of images.

- a.) Make ONE aesthetically pleasing PPT slide or jupyter notebook from all images acquired in your laboratory course section.
- b.) Calculate approximate magnification of objective lens (lens gap is 2.2 mm and sample sits in the middle)

5 IMAGE ANALYSIS SOFTWARE FOR SEM AND TEM

- **Digital Micrograph:** You can get a free offline license from here:
<https://www.gatan.com/resources/scripting/demo/>
- **ImageJ:** a free imaging software with a lot of plugins download from here:
<https://rsbweb.nih.gov/ij/>
- **jupyter notebook:** use provided jupyter notebook
<https://github.com/gduscher/MSE672-Introduction-to-TEM>